



## Bair Hugger® Therapy Product Specifications

### BAIR HUGGER BLANKETS

Bair Hugger therapy offers 25 styles of blankets to meet both your everyday and specialized patient warming needs—from pediatric to geriatric, from brief outpatient procedures to long complex procedures.


#### All Bair Hugger blanket features:

- Flexible, easy to use designs
- Uniform perforation pattern across the blanket surface ensures even convective warming
- Resealable hose ports (where applicable)
- Soft, radiolucent, latex-free materials

### INTRAOPERATIVE BLANKETS: UPPER, LOWER, FULL & TORSO


#### Upper Body Blanket - Model 522

<b>Size:</b> 74 x 24 in (188 x 61 cm)	<b>Weight:</b> 3.5 oz (99 g)	<b>Drape:</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 2
--	---------------------------------	--	-------------------------



#### XL Upper Body Blanket - Model 523

<b>Size:</b> 84 x 36 in (213 x 91 cm)	<b>Weight:</b> 5.4 oz (153 g)	<b>Drape:</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 2
--	----------------------------------	--	-------------------------



#### Lower Body Blanket - Model 525

<b>Size:</b> 60 x 36 in (152 x 91 cm)	<b>Weight:</b> 4.0 oz (113 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
--	----------------------------------	----------------------	-------------------------



#### Torso Blanket - Model 540

<b>Size:</b> 42 x 36 in (107 x 91 cm)	<b>Weight:</b> 3.2 oz (90 g)	<b>Drape:</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 1
--	---------------------------------	--	-------------------------



#### Dual Port Torso Blanket - Model 542

<b>Size:</b> 42 x 36 in (107 x 91 cm)	<b>Weight:</b> 3.2 oz (90 g)	<b>Drape:</b> 24 x 16 in (61 x 41 cm)	<b>Hose Ports:</b> 2
--	---------------------------------	--	-------------------------





## Bair Hugger® Therapy Product Specifications

### UNDERBODY SERIES BLANKETS

Prewarm your patients - As little as 15 minutes of pre-induction forced-air warming can reduce the severity of re-distribution temperature drop. Now you can warm your patients from the start.

#### The Bair Hugger Underbody Series blanket features:

- Advanced blanket designs warm the patient while providing full, unrestricted patient access
- Patented fluid outlets minimize pooling of fluids on the surface of the blanket
- Pass-through slit design allows for flexible patient positioning and accommodates the use of a drawsheet (where applicable)
- Soft, radiolucent, latex-free materials
- Resealable hose ports at either end of the blanket provide options for hose placement (where applicable)

#### Pediatric Underbody Blanket - Model 555

Size:	Weight:	Drape (2):	Hose Ports:
36 x 33 in (91 x 84 cm)	3.0 oz (85 g)	24 x 24 in (61 x 61 cm)	2



#### Large Pediatric Underbody Blanket - Model 550

Size:	Weight:	Drape (2):	Hose Ports:
60 x 32 in (152 x 81 cm)	4.8 oz (136 g)	24 x 24 in (61 x 61 cm)	2



#### Adult Underbody Blanket - Model 545

Size:	Weight:	Drape:	Hose Ports:
74 x 36 in (188 x 91 cm)	5.0 oz (142 g)	n/a	1




#### Spinal Underbody Blanket - Model 575

Size:	Weight:	Drape:	Hose Ports:
91 x 36 in (213 x 91 cm)	8 oz (218 g)	24 x 48 in (61 x 122 cm)	1



#### Lithotomy Underbody Blanket - Model 585

Size:	Weight:	Drape (1):	Hose Ports:
64 x 36 in (160 x 91 cm)	5.4 oz (153 g)	24 x 24 in (61 x 61 cm)	1



#### Full Access Underbody Blanket - Model 635

Size:	Weight:	Drape (1):	Hose Ports:
84 x 36 in (221 x 91 cm)	7.0 oz (198 g)	24 x 48 in (61 x 122 cm)	2, presealed



#### Sterile Full Access Underbody Blanket - Model 637

Size:	Weight:	Drape (1):	Hose Ports:
84 x 36 in (221 x 91 cm)	7.0 oz (198 g)	24 x 48 in (61 x 122 cm)	2, presealed





## Bair Hugger® Therapy Product Specifications

### PEDIATRIC BLANKETS

Bair Hugger therapy provides effective temperature management for the premature infant, neonate or child in any clinical setting. There is no longer a need to adapt products and equipment for pediatric use. Our pediatric warming solutions are designed specifically with younger patients in mind.

#### The Bair Hugger Pediatric blankets feature:

- Flexible, easy to use designs
- Uniform perforation pattern across the blanket surface ensures even convective warming
- Soft, radiolucent, latex-free materials

#### Small Lower Body Blanket - Model 537

<b>Size:</b> 35 x 24 in (89 x 61 cm)	<b>Weight:</b> 1.5 oz (43 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
---	---------------------------------	----------------------	-------------------------



#### Pediatric Long Blanket - Model 530

<b>Size:</b> 74 x 8.5 in (188 x 21 cm)	<b>Weight:</b> 10 oz (28 g)	<b>Drape (2):</b> 24 x 18 in (61 x 46 cm)	<b>Hose Ports:</b> 1
---	--------------------------------	--	-------------------------



#### Pediatric Underbody Blanket - Model 555

<b>Size:</b> 36 x 33 in (91 x 84 cm)	<b>Weight:</b> 3.0 oz (85 g)	<b>Drape (2):</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 2
---	---------------------------------	--	-------------------------



#### Large Pediatric Underbody Blanket - Model 550

<b>Size:</b> 60 x 32 in (152 x 81 cm)	<b>Weight:</b> 4.8 oz (136 g)	<b>Drape (2):</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 2
--	----------------------------------	--	-------------------------



#### Postoperative Pediatric Full Body Blanket - Model 310

<b>Size:</b> 60 x 36 in (152 x 91 cm)	<b>Weight:</b> 4.8 oz (134 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
--	----------------------------------	----------------------	-------------------------





## Bair Hugger® Therapy Product Specifications

### SPECIALTY & CARDIAC BLANKETS

Challenging clinical scenarios require unique solutions. Bair Hugger therapy offers seven unique styles of specialty and cardiac blankets to help you reach your warming goals.

#### The Bair Hugger Specialty & Cardiac blankets feature:

- Flexible, easy to use designs
- Uniform perforation pattern across the blanket surface ensures even convective warming
- Soft, radiolucent, latex-free materials

#### Preoperative & Outpatient Care Blanket - Model 110

<b>Size:</b>	<b>Weight:</b>	<b>Booties:</b>	<b>Hose Ports:</b>
84 x 36 in (213 x 91 cm)	5.6 oz (157 g)	1 pair	1



#### Preoperative & Outpatient Care Blanket - Model 111

<b>Size:</b>	<b>Weight:</b>	<b>Booties:</b>	<b>Hose Ports:</b>
84 x 36 in (213 x 91 cm)	4.8 oz (134 g)	n/a	1



#### Cath Lab Blanket - Model 560

<b>Size:</b>	<b>Weight:</b>	<b>Drape:</b>	<b>Hose Ports:</b>
132 x 9 in (335 x 23 cm)	3.0 oz (85 g)	n/a	1



#### Surgical Access Blanket - Model 570

<b>Size:</b>	<b>Weight:</b>	<b>Drape:</b>	<b>Hose Ports:</b>
84 x 36 in (213 x 91 cm)	6.4 oz (181 g)	24 x 16 in (61 x 41 cm)	2



#### Full Body Surgical Blanket - Model 610

<b>Size:</b>	<b>Weight:</b>	<b>Drape:</b>	<b>Hose Ports:</b>
72 x 36 in (183 x 91 cm)	4.8 oz (134 g)	n/a	1



#### Sterile Cardiac Blanket - Model 630

<b>Size:</b>	<b>Weight:</b>	<b>Drape:</b>	<b>Hose Ports:</b>
60 x 36 in (152 x 91 cm)	8.0 oz (227 g)	n/a	1



#### Sterile Cardiac Access Blanket - Model 645

<b>Size:</b>	<b>Weight:</b>	<b>Drape:</b>	<b>Hose Ports:</b>
108 x 9 in (274 x 23 cm)	4.0 oz (113 g)	n/a	1





## Bair Hugger® Therapy Product Specifications

### PACU BLANKETS

Bair Hugger therapy is as much at home postoperatively as it is in the operating room. A wide range of blanket styles gives you optimum flexibility to treat virtually any patient.

#### The Bair Hugger Postoperative Care blanket features:

- Patented blanket designs that provide even heat distribution
- Durable and flexible materials that provide the ability to warm patients in various positions
- Non-inflated foot drapes guard against thermal injury to the heat-sensitive feet and lower legs

#### Full Body Blanket - Model 300

<b>Size:</b> 84 x 36 in (213 x 91 cm)	<b>Weight:</b> 5.1 oz (145 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
--	----------------------------------	----------------------	-------------------------



#### Chest Access Blanket - Model 305

<b>Size:</b> 72 x 36 in (183 x 91 cm)	<b>Weight:</b> 4.7 oz (133 g)	<b>Drape (2):</b> 24 x 24 in (61 x 61 cm)	<b>Hose Ports:</b> 1
--	----------------------------------	--	-------------------------



#### Postoperative Pediatric Full Body Blanket - Model 310

<b>Size:</b> 60 x 36 in (152 x 91 cm)	<b>Weight:</b> 4.8 oz (134 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
--	----------------------------------	----------------------	-------------------------



#### Multi-Access Blanket - Model 315

<b>Size:</b> 84 x 36 in (213 x 91 cm)	<b>Weight:</b> 5.3 oz (150 g)	<b>Drape:</b> n/a	<b>Hose Ports:</b> 1
--	----------------------------------	----------------------	-------------------------





## Bair Hugger® Therapy Product Specifications

### TEMPERATURE MANAGEMENT UNITS

Bair Hugger® therapy temperature management units have become the standard of excellence in forced-air warming. They are small, lightweight and easy to use. All Bair Hugger temperature management units are designed to provide clinically effective patient warming at safe air delivery temperatures.

#### Temperature Management Unit - Model 505

**Dimensions:** 13" h x 10" w x 11" d (33 x 25 x 28 cm)

**Weight:** 13.6 lb (6.2 kg)

**Filter:** High-efficiency, 0.2 µm filter

**Operating Temperature:**

High: 43° ± 3°C (109.4° ± 5.4°F)

Med: 38° ± 3°C (100.4° ± 5.4°F)

Low: 32° ± 3°C (89.6° ± 5.4°F)

**Device Rating:**

110-120 VAC, 60 Hz, 9.5 Amperes

220-240 VAC, 50 Hz, 4.5 Amperes

100 VAC, 50/60 Hz, 10 Amperes

**Certifications:** UL 544, CSA C22.2 No. 125, IEC 60601-1-1, IEC 60601-1-2, EN55014, AS 3200.1990

**Leakage Current:** Meets regulatory standards for leakage current



#### Temperature Management Unit - Model 750

**Dimensions:** 12" h x 10" w x 13.5" d (30 x 25 x 34 cm)

**Weight:** 15.5 lb (7 kg)

**Filter:** High-efficiency, 0.2 µm filter

**Operating Temperature:**

High: 43° ± 1.5°C (109.4° ± 2.7°F)

Med: 38° ± 1.5°C (100.4° ± 2.7°F)

Low: 32° ± 1.5°C (89.6° ± 2.7°F)

**Device Rating:**

110-120 VAC, 50/60 Hz, 11.7 Amperes

220-240 VAC, 50/60 Hz, 7.2 Amperes

100 VAC, 50/60 Hz, 14.5 Amperes

**Certifications:**

EN 60601-1, UL 2601-1, CSA C22.2 No. 601.1, EN 60601-1-2

**Leakage Current:** Meets regulatory standards for leakage current

